

**Title:**

Reverse Transcription (RT) for First-strand cDNA synthesis

Date last updated:

20<sup>th</sup> June 2024

Total Pages:

2

Written by:

Jeanie Wu

**Reverse Transcription (RT) for First-strand cDNA synthesis****I. PURPOSE**

To reverse transcribe total RNA and synthesize first-strand cDNA for downstream Next-Generation Sequencing library preparation.

**II. MATERIALS & EQUIPMENTS**

Sample type

- Extracted arboviral RNA

Reagents

- SuperScript® III First-Strand Synthesis System for RT (Invitrogen Cat# 18080051)
- Nuclease – free ultrapure water

Consumables

- 0.2ml PCR tubes

Equipment

- Thermal cycler

**III. PROCEDURE**

Note: Use 8ul (maximum amount) of extracted RNA and 1ul of random hexamers for this protocol with Dengue and Zika virus.

1. Prepare RNA mixture:

## a) Prepare RNA mixture with the following components in a 0.2 ml PCR tube

Component	Amount/ sample
Extracted Sample RNA	8 uL
Primers - 50 ng/μL random hexamers	1 uL
10 mM dNTP mix	1 uL
Nuclease-free water *	0 ul
<b>Total volume</b>	<b>10 ul</b>

\* top up to 10 ul water if sample is less than 8 ul

**Title:**

Reverse Transcription (RT) for First-strand cDNA synthesis

Date last updated:  
20<sup>th</sup> June 2024Total Pages:  
2Written by:  
Jeanie Wu

- b) Place tube in thermal cycler at 65°C for 5 mins
- c) Then incubate on ice for 1 min and quick spin down.

- Prepare the cDNA synthesis mix.

- a) In a separate tube, add each component in the following indicated order:

Component	1 Rxn	10 Rxns
10X RT buffer	2 uL	20 uL
25 mM MgCl <sub>2</sub>	4 uL	40 uL
0.1 M DTT	2 uL	20 uL
RNaseOUT™ (40 U/μL)	1 uL	10 uL
SuperScript® III RT (200 U/μL)	1 uL	10 uL
<b>Total volume</b>	<b>10 uL</b>	<b>100 uL</b>

- b) Mix well by gently pipetting up and down. Quick spin down.

- Pipette 10ul of cDNA synthesis mix to each tube with sample RNA mixture. Gently pipet to mix and quick spin.
- Place PCR tubes in thermocycler and incubate as follows:
  - a. Random hexamer primed: 25°C for 10 mins, followed by 50°C for 50 mins
- Inactivate the reaction by heating to 85°C for 5 mins. Then chill on ice.

***Note: The cDNA can now be used as a template for amplification in PCR. However, amplification of some PCR targets (those > 1 kb) may require the removal of RNA complementary to the cDNA.***

- To remove RNA complementary to cDNA, add 1 ul (2units) of *E. coli* RNase H.
- Incubate at 37°C for 20 mins.
- cDNA synthesis reaction can be stored at -20°C or used for PCR immediately.