

Introduction to primer design concepts and tools

“Arboviral Genomic Surveillance in a Climate evolving world” Next Generation Sequencing (NGS) workshop


Duke-NUS Medical School

Singapore

Presented by: Dr. Chanditha Hapuarachchi

03 Jul 2024

Next Generation Sequencing for genomic surveillance



**Target
enrichment**

- **PCR**
- Isolation

Long oligo
probes


Background
depletion

Amplicon sequencing

Amplicon sequencing is a type of **targeted** next generation sequencing that uses **multiplex PCR** to create DNA sequences called amplicons

Advantages



- Fast turnaround time
- Low sequencing costs
- Low DNA input requirements – up to 50-100 genome copies per reaction
- Supports multiplexing of multiple samples in a single sequencing reaction
- Amenable to a wide range of experimental designs, such as **complete genome sequencing, allelic variant identification, marker gene sequencing for microbiome studies** etc.



A successful PCR relies on well-designed primers

What makes a good primer?

High priming efficiency	High specificity
No homo/hetero dimers, esp. at 3' end (prevent self-extension)	Long enough and reasonably high melting temperature (T_m)
No hairpin stems esp. at 3' end (prevent self-priming)	Uniqueness of the sequence
No repetitive sequences & mismatches esp. at 3' end	Moderately stable 3' end



Characteristics of a “good” primer



Requirement	Feature
High priming efficiency	Approximate T_m values ($\pm 5^\circ\text{C}$ between F & R)
	Equal proportions of A/Ts & G/Cs (50% $\pm 5\text{-}10\%$)
	Includes G/Cs at the 3' end
	Lack of long repeats & runs— no more than 4 repeats & 3 bp runs
	Lack of duplexes and structures
High specificity	Length 18-22 bp
	Unique sequence esp. at 3' end but with moderate stability (ideally no more than 2 Gs or Cs in last 5 bases)

Other considerations

- Avoid regions with high secondary structures (hairpins) & GC islands
- Try to have the lowest level of degeneracy, if degenerate bases are needed
- **For difficult to amplify regions**
 - Select primers with high T_m and GC content (than usual)
 - Add solvents/additives such as DMSO
 - Increase primer concentration, but not more than 1 μmol (to avoid nonspecific primer-primer amplifications)
- **Designing probes**
 - Do not put Gs at the 5' end
 - Keep the size of PCR product or amplicon to 50 - 150 bp. Shorter amplicons work most efficiently
 - Probe should ideally be close to one of the primers without overlapping

Post-design analysis

OligoAnalyzer

Sequence

5' MOD INTERNAL 3' MOD MIXED BASES

TACAACATGATGGGGAARAGAGARAA

Bases 26

CLEAR SEQUENCE

Try the new batch mode here

Parameter sets

SpecSheet (Default)

Target type

DNA

Oligo Conc

0.25 μ M

Na⁺ Conc

50 mM

Mg⁺⁺ Conc

1.50 mM

dNTPs Conc

200 mM

Instructions | Definitions | Feedback

ANALYZE

HAIRPIN

SELF-DIMER

HETERO-DIMER

NCBI BLAST

TM MISMATCH

ADD TO ORDER


Results

RESUSPENSION DILUTION

SEQUENCE	5'- TAC AAC ATG ATG GGG AAR AGA GAR AA -3'		
COMPLEMENT	5'- TTY TCT CTY TTC CCC ATC ATG TTG TA -3'		
LENGTH	26		
GC CONTENT	38.5 %		
MELT TEMP RANGE	MIN	MEAN	MAX
	63 °C	63.6 °C	64.7 °C

Structures

structure	Image	ΔG (kcal.mole ⁻¹)	T _m (°C)	ΔH (kcal.mole ⁻¹)	ΔS (cal.K ⁻¹ .mole ⁻¹)	Output
1		-0.04	25.5	-21.9	-73.33	<div>Ct</div> <div>Det</div>
2		0.02	24.6	-14.7	-49.38	<div>Ct</div> <div>Det</div>
3		0.41	17.7	-16.5	-56.73	<div>Ct</div> <div>Det</div>



Homo-Dimer Analysis

The delta G is calculated by taking into account the longest stretch of complementary bases. These pairs of complementary bases are represented by a solid line. Dotted lines represent additional complementary bases for that dimer structure, but their presence does not impact calculated delta G values. Actual delta G values may vary based on presence of additional complementary bases. The Maximum Delta G value refers to the free energy of the oligo sequence binding to its perfect complement.

Dimer Sequence:
5'- TACAACATGATGGGGAARAGAGARAA -3'
Maximum Delta G: -46.01 kcal/mole

Delta G: -5.38 kcal/mole Base Pairs: 4
5' TACAACATGATGGGGAARAGAGARAA
:: |||| ::
3' AARAGAGARAAGGGGTAGTACAACAT

Hetero-Dimer Analysis


The delta G is calculated by taking into account the longest stretch of complementary bases. These pairs of complementary bases are represented by a solid line. Dotted lines represent additional complementary bases for that dimer structure, but their presence does not impact calculated delta G values. Actual delta G values may vary based on presence of additional complementary bases. The Maximum Delta G value refers to the free energy of the oligo sequence binding to its perfect complement.

Primary Sequence: 5'- TACAACATGATGGGGAARAGAGARAA -3'
Secondary Sequence: 5'- GTGTCCCAGCCGCGTGCATCGC -3'
Maximum Delta G: -50.44 kcal/mole

Delta G: -8.09 kcal/mole Base Pairs: 4
5' TACAACATGATGGGGAARAGAGARAA
::: : |||| :
3' CGCTACTGTGCGCCGACCCCTGTG

Primer designing for amplicon sequencing

- 'Primal Scheme' primer design tool (<https://primalscheme.com/>)
- Multiplex primer sets that are divided into two pools (pool 1 & 2)
- User-defined fragment length (default 400 bp; 200-2000 bp)
- Overlapping fragments – tiling along the genome with overlapping regions (default 75 bp)
- High annealing temperature that is high enough to perform 2-step PCR (65-68°C)



The screenshot shows the 'primalscheme' web interface. The header includes the logo 'primalscheme' and the tagline 'primer panels for multiplex PCR'. A 'Design a new scheme' form is displayed with a 'Reset defaults' button. The form contains a 'FASTA file...' upload section, 'Options' with 'High-GC mode' and 'Pinned' toggles, and an 'Amplicon size' input set to 400. A 'Scheme name' input is set to 'e.g., nCov-400'. A 'Design scheme' button is at the bottom, with a 'View demo inputs' link below it. A small blue and yellow starburst icon with '00' is in the top right corner of the interface.

Quick J et al. Multiplex PCR method for MinION and Illumina sequencing of Zika and other virus genomes directly from clinical samples. Nat Protoc. 2017 Jun;12(6):1261-1276.

Example panels

1. DENV (Nathan Grubaugh)


https://www.protocols.io/view/denguese-a-pan-serotype-whole-genome-amplicon-seq-kqdq39xxeg25/v2?version_warning=no

2. ZIKV & CHIKV

Quick J et al. Multiplex PCR method for MinION and Illumina sequencing of Zika and other virus genomes directly from clinical samples. Nat Protoc. 2017 Jun;12(6):1261-1276.

Primalscheme: Primer designing workflow

Design requirement - research question, target region etc



Primalscheme: Primer designing output files


BED file – Provides primer coordinates

```
AY618991(country      39      61      DENV-4_Revised_1_LEFT  1      +
AY618991(country      518     541     DENV-4_Revised_1_RIGHT 1      -
AY618991(country      444     466     DENV-4_Revised_2_LEFT  2      +
AY618991(country      904     926     DENV-4_Revised_2_RIGHT 2      -
AY618991(country      831     852     DENV-4_Revised_3_LEFT  1      +
AY618991(country     1312    1334    DENV-4_Revised_3_RIGHT 1      -
```

TSV file – Provides primer sequences

name	pool	seq	size	%gc	tm (use 65)
DENV-4_Revised_1_LEFT	1	ATGCTGAAACGCGAGAGAAACC	22	50	61.93
DENV-4_Revised_1_RIGHT	1	ACCAGCAATCAATGTCTTCAGGT	23	43.48	60.76
DENV-4_Revised_2_LEFT	2	CTTATTGCCATGGACCTGGGTG	22	54.55	61.46
DENV-4_Revised_2_RIGHT	2	CATCCTCCATGTTCTAGCACCA	22	50	60.01
DENV-4_Revised_3_LEFT	1	TACGGAATGCGATGCGTAGGA	21	52.38	61.79
DENV-4_Revised_3_RIGHT	1	GGGGTTATCGTAGCTGTCACTC	22	54.55	60.4


PDF file – tiling plot




Amplicon sequencing: Primer designing considerations

- Selection of sequences: to represent the true diversity of the pathogen/genomic region of interest
- No. of panels: based on the level of genetic divergence
- Quality of nucleic acids: Short amplicons for degraded samples

High variant diversity within a single clade of DENV-2




Dengue virus diversity cascade: Ideally needs panels to represent each genotype



Amplicon sequencing: Primer designing troubleshooting


Primer panels do not always generate ideal results!!

Drop-out regions




Design limitations

500 bp design





2000 bp design



Amplicon sequencing: Primer designing troubleshooting

Improving the primers in drop-out regions



- Map the primers to determine the likely reason of their failure
 - Modify the primer if possible, while maintaining the same T_m
 - Choose a replacement primer from another panel
 - Redesign the panel using a different reference dataset
- If fails, fill-in gaps using Sanger sequencing

Our Environment

Safeguard • Nurture • Cherish